

Bringing Together the Best in Laparoscopic Surgical Instrumentation

Handpieces 5mm Limited Use Monopolar



- 3904 ReNew IV Handpiece 34cm
- 3905 ReNew IV Handpiece 42 cm
- 3908 ReNew IV Handpiece 25cm
- 3924 ReNew RS BARIATRIC Handpiece 34cm
- 3925 ReNew RS BARIATRIC Handpiece 42cm




- 3914 ReNew IV Ratcheted Handpiece 34cm
- 3915 ReNew IV Ratcheted Handpiece 42 cm
- 3918 ReNew IV Ratcheted Handpiece 25cm
- 3934 ReNew RS BARIATRIC Ratcheted Handpiece 34cm
- 3935 ReNew RS BARIATRIC Ratcheted Handpiece 42 cm



Tips

EverSharp® Disposable Scissors



Disposable ControlTip® Graspers and Dissectors



Disposable ControlTip® Super Atrau Graspers



Limited Use ControlTip® Graspers and Dissectors



Limited Use EverSharp® Scissors



© 2005 ControlTip®, EverSharp®, and ReNew™ are trademarks of Microline Surgical, Inc. All rights reserved. P/N 09-39-00064-00, 12/09

Confidence simply delivered

Tip connection Tips!

- 1 Make sure the protective white O-ring is fully intact
- 2 Microline's patented co-axial threading ensures easy and safe connection
- 3 Close handle and tip then connect clockwise ...
- 4 ... ensure the white O-ring is fully invisible, for full HF safety



Authorized EU Representative:
Medical Device Safety Service GmbH
Schiffgraben 41 - 30175 Hannover - Germany
www.microlinesurgical.com

